Monte Carlo Simulations of a RENS Operating ToF Neutron Spectrometer

A. Cannuli¹, M.T. Caccamo¹, S.Magazú¹

E-mail: acannuli@unime.it

Abstract

In the present work the results of simulations performed by means of the McStas Monte Carlo simulator software package are presented. In particular simulation results obtained with modified instrumental parameters, i.e. energy resolution and wave-vector range, through McStas interface allowed to evaluate the IN5 instrument behavior when operating following the RENS protocol.

Keywords: Resolution Elastic Neutron Scattering (RENS), McStas MonteCarlo simulations, European Spallation Source (ESS).

Introduction

The aim of this work is to simulate the Resolution Elastic Neutron Scattering (RENS) approach [1]-[2], which is based on the acquisition of the measured elastic scattering law as a function of the instrumental energy resolution, for the IN5 instrument, a general purpose direct geometry Time of Flight (ToF) spectrometer of the Institut Laue- Langevin (ILL) of Grenoble, France. In the RENS approach the measured elastic scattering law as a function of the logarithm of the instrumental energy resolution behaves like an increasing sigmoid curve whose inflection point occurs when this instrumental resolution time matches the system relaxation time. Simulation results performed by means of the McStas Monte Carlo simulator software package to obtain the behavior of the IN5 instrument operating following the RENS protocol are presented. At first, a brief description of IN5 instrument is dealt; then, the simulations performed assuming a wavelength of 5, 10 and 15 \mathring{A} at different speeds of rotation of the chopper respectively: 2000, 5000, 8000, 11000, 14000 and 17000 rpm. In particular, the first simulation with dynamical structure factor for Vanadium, obtained by keeping constant the thickness of the sample 1 mm and the incoming wavelength spread $\Delta\lambda$ 3 Å, the wavelength of 5 \mathring{A} and the chopper speed of 11000 rpm, is reported. The final result of all simulations is the IN5 two-dimensional wavelength-choppers speed plane that confirm the validity of the RENS approach which is based on the acquisition of the measured elastic scattering as a function of the instrumental energy resolution.

Methods: McStas software package and IN5 spectrometer

McStas is a Monte Carlo simulation software, programmed in C and assisted by Phyton-code for the interface (Graphic User Interface (GUI)). Born from the collaboration between DTU Physics (RISO DTU and RISO Natnl. Lab), University of Copenhagen, Paul Scherrer Institute and ILL (Institut Laue-Langevin) in October 10, 2001; the current version is 2.1, released on 12-09-2014. The Monte Carlo method is used for the management of neutron interaction: the system generates a random number which is compared with a probability and, depending on the number extract, the particle assumes a given behavior [3]-[4]. It possible visualize 3D graphic of the desired simulated instrument. It allows to simulate what happens changing of parameters (e.g. geometry, chopper velocity, etc) and to evaluate the best conditions of use for existing instruments and to design new spectrometers. The present version supports both continuous and pulsed source instruments, and includes a library of standard components with in total around 130 components. IN5 is a general purpose direct geometry Time of Flight (ToF) spectrometer. It is used to study low-energy transfer processes as a function of momentum transfer. Typically this instrument is used for measurements

¹ Dipartimento di Scienze Matematiche e Informatiche, Scienze Fisiche e Scienze della Terra, viale F. Stagno d'Alcontres 31. Messina. Italy

in the small energy and momentum transfer region with values of the order of $\delta E/E = 1$




Figure 1: 3D simulation of IN5 instrument

Primary spectrometer: the chopper system provides an extremely well defined resolution function which remains nearly Gaussian (counter-rotating disks) down to 1/10.000 of its height. Secondary spectrometer: it consists in a large cylindrical array of pixelated position sensitive detector (PSD) mounted inside a vacuum time-of-flight chamber 4m from the sample. PSDs cover a large gapless detection surface (30m²) with a spatial resolution comparable to triple-axis instruments. With the new detector, the instrument has a large solid-angle for powder and amorphous systems experiments and is fully equipped for single-crystal spectroscopy. Sample area: A large sample chamber (max diameter = 80cm, depth ~ 130 cm) allows hosting the broad range of IN5 sample environments but also user-designed equipment (contact instrument responsible before). The non-magnetic material surrounding the sample allows the use of either high magnetic field (limited to 2.5T at the moment it is foreseen to purchase a specific magnet for higher fields - equal or above 10T) or neutron polarization devices (still to come).




Figure 2: 2D simulation of IN5 instrument

IN5 is constituted by a series of guides, chopper and collimators, which accompany the neutron beam from the source to the sample. The neutrons diffuse along a series of optical guides, which avoids the dispersion of the beam, their path is 17 meters. The chopper are 7: the chopper 0, 1, 2, 4, 5 and 6 all rotate at the same speed, the 3 has a speed proportional to that of the other chopper, usually half and this allows us to increase the interval between one

pulse and another. The chopper 0, 1, and 2 have an opening angle of 9°, on 3 and 4 to 9.5°, 5 and 6 have angular aperture of 3.25°. A small opening has the purpose of monochromate the neutron beam. The guides have a total length of about 15 meters and are intended to focus the beam, thus to prevent the dispersion. The collimator has the purpose of making the neutron beam parallel, we IN5 of about 1 meter of collimator placed before the sample.

McStas IN5 experimental results

We have performed the following simulations where some parameters were kept constant whereas r rotation speed within the limits of operation. Such limit parameters are reported in Table IN5 of the preceding paragraph. The simulations were performed assuming a wavelength of 5, 10 and 15 \mathring{A} at different speeds of rotation of the chopper respectively: 2000, 5000, 8000, 11000, 14000 and 17000 rpm. More precisely in the first simulation the instrument parameters are: i) the thickness of the sample; having chosen a cylindrical can, taking into account both the sides of cylinder, the total thickness was 1 mm; ii) the incoming wavelength spread $\Delta\lambda$, which was kept to the value of 3 \mathring{A} ; iii) the simulation was performed assuming a wavelength of 5 \mathring{A} ; iv) the chopper speed is 11000rpm. The chosen set-up parameter values of such a simulation are reported in Figure 3.

INSTRUMENT'S SETTING					
lambda	5 Å	dlambda	зъ	speed	11000rpm
ratio	0.5	housing	Fe.laz	coh	V.laz
inc	null	thickness	1 mm	height	2.5 mm
radius	5 mm		order	multiple scattering	
Neutron count			50000000		

Figure 3: The table synthesizes the employed fitting parameters of the performed simulation. In particular, the constant instrument parameters are: i) the thickness of the sample: 1 mm; ii) the incoming wavelength spread $\Delta\lambda$: 3 Å.

The obtained output of the performed simulation results a Gaussian like curve whose width has been extracted through a fitting procedure obtained by means of the commercially available software ORI-GIN 9.0. The obtained Full Width at Half Maximum (FWHM) is: ΔE 79.89 μeV . In the following picture the obtained result is reported withing the two-dimensional plot wavelength-choppers speed. As it can be seen the agreement is very satisfactionable.




Figure 4: Simulated dynamical structure factor for Vanadium as obtained by keeping constant the thickness of the sample: 1 mm and the incoming wavelength spread $\Delta\lambda$: 3 Å. The incoming wavelength is: 5 Å; the chopper speed is 11000 rpm; the neutron counts are fixed to 500.000.000; the speed ratio between chopper 3 and chopper 1 is 0.5. The obtained Gaussian like curve shows Full Width at Half Maximum (FWHM) of ΔE 79.89 μeV

The results of all simulations are reported below, in the IN5 two-dimensional wavelength-choppers speed plane. The points in yellow represents the values of the obtained resolutions at the following experimental conditions: i) incoming wavelength =5, 10 and 15 Å; chopper speed =2000, 5000, 8000, 11000, 14000 and 17000 rpm. The obtained FWHMs are:

$\Delta E = 432.95 \ \mu eV;$	$\Delta E = 10.42 \ \mu eV;$
$\Delta E = 174.83 \ \mu eV;$	$\Delta E =~8.43~\mu eV;$
$\Delta E = 110.38 \ \mu eV;$	$\Delta E = 7.16 \mu eV;$
$\Delta E = 79.89 \mu eV;$	$\Delta E = 16.46 \ \mu eV;$
$\Delta E = 62.04 \mu eV;$	$\Delta E = 6.74 \ \mu eV;$
$\Delta E = 50.74 \mu eV;$	$\Delta E = 4.27 \mu eV;$
$\Delta E = 55.82 \mu eV;$	$\Delta E = 3.13 \ \mu eV;$
$\Delta E = 22.16 \mu eV;$	$\Delta E = 2.05 \ \mu eV;$
$\Delta E = 14.02 \mu eV;$	$\Delta E = 2.06 \mu eV$.

Figure 5: Obtained Full Width at Half Maximum (FWHM) values




Figure 6: The IN5 two-dimensional wavelength-choppers speed plane.

Conclusions

In this paper McStas Monte Carlo simulations addressed to simulate the IN5 instrument behavior has been performed. A varying instrumental energy resolution, and a varying energy window as well, allow to characterize the RENS approach. Overall, this work confirms the validity of the RENS approach which is based on the acquisition of the measured elastic scattering as a function of the instrumental energy resolution; in particular, following it, the measured elastic scattering law as a function of the logarithm of the instrumental energy resolution behaves like an increasing sigmoid curve whose inflection point occurs when this the instrumental resolution time matches the system relaxation time.

References

- [1] Magazú, S. et al.: Elastic Incoherent Neutron Scattering Operating by Varying Instrumental Energy Resolution: Principle, Simulations and Experiments of the Resolution Elastic Neutron Scattering (RENS). Rev. Sci. Instr.., 82, 2011, p. 9.
- [2] Magazú, S., Migliardo, F., Caccamo, M.T.: Upgrading of Resolution Elastic Neutron Scattering (RENS), Advances in Materials Science and Engineering., 2013, and references therein.
- [3] Cannuli, A.: Resolution Elastic Neutron Scattering spectrometer operating through a step-by-step rotating monochromato Thesis.

[4] Caccamo, M.T., Magazú, S.: McStas simulation on IN5 spectrometer following RENS approach.

volume: Activity Report 2014, ISSN 2038-5889, ed. L. Torrisi (2014).